

Dräger



DEUTSCHER TEXT:
BITTE UMDREHEN

OPERATING MANUAL

Romulus® 800 Anaesthetic Apparatus

A 14

FROM Dräger.

Romulus® 800/800 V/800 M/800 MV

OPERATING INSTRUCTIONS

Important Notice

For correct and effective use of the device, and to avoid hazards, we would point out the following:

- 1 Any use of the device requires precise knowledge and observation of these operating instructions.
- 2 The device is intended only for the purposes specified in the Operating Manual or for purposes confirmed in writing by Drägerwerk AG.
- 3 The device should be inspected by experts at regular time intervals. An official report of the inspections should be drawn up.
- 4 Only original Dräger spare parts should be used for maintenance and repairs. Repairs and maintenance, and the replacement of spare parts should only be carried out by experts.
- 5 We recommend having inspections and repair work carried out by the Technical Customer Service of your Dräger Branch or Agent.

Regular inspection is best ensured by entering into an Inspection Service Contract with the Technical Customer Service of your Dräger Branch or Agent.

- 6 Responsibility for the reliable function of the device passes to the owner or operator in all cases where the device has been inexpertly maintained or repaired by persons not employed by the Dräger Organization or where it has been used in a manner which does not conform to the normal conditions of use.
- 7 For reasons of safety, pressure reducers should be overhauled at least every 6 years.

We would also point out that the national recommendations, regulations and laws governing the use of technical equipment should be observed.

DRÄGERWERK AG LÜBECK

Contents

	Page		Page		Page
Important Notice	2	5 Testing Readiness for Operation	13	6.7 Bag for manual ventilation .	16
1 Intended Use	3	5.1 Gas supply	13	7 Shut-Down Actions	16
1.1 Applications	3	5.2 Gas deficiency signal, blender function	13	8 Care	16
1.2 Notes on safety	3	5.3 O₂ flush (bypass)	14	8.1 Stripping down	16
2 Design and Function	6	5.4 Circle system	14	8.2 Hygiene	16
3 Alarm and Safety Devices	8	5.5 Ventilog	14	8.3 Assembly	17
3.1 O₂ deficiency signal, N₂O cutoff	8	5.6 Secretion aspirator	14	9 Servicing	17
3.2 Other safety features	9	5.7 Microbe filter	14	10 Technical Data	17
4 Initial Preparation	9	5.8 Vapor 19.1	14	11 Description of Special Versions Romulus 800 »Air« / 800 V »Air«	19
4.1 Gas supply	9	5.9 Anaesthetic gas removal	14	12 Check List	22
4.2 Ventilog	10	5.10 Bag for manual ventilation .	14	13 Order List	24
4.3 Circle system	10	6 Operational Use	15	14 Parts List	26
4.4 Secretion aspirator	11	6.1 Possible ventilation modes .	15		
4.5 Microbe filter	11	6.2 Secretion aspiration	15		
4.6 Vapor 19.1	12	6.3 Vapor 19.1	15		
4.7 Anaesthetic gas removal	12	6.4 O₂ flush (bypass)	15		
4.8 Potential equalization	12	6.5 O₂ deficiency signal	15		
4.9 Bag for manual ventilation .	12	6.6 N₂O deficiency	16		

Applicable model designation:
 (see plate on front of device)

Fabrication No.:
 (next to model designation)

Explanation of model designations

Romulus Model	Gas blending by means of	Ventilator Ventilog®	Remarks
800	Flowmeter unit ¹⁾	Attachment possible	
800 V	Flowmeter unit ¹⁾	Integrated	
800 M	Gas blender ²⁾	Attachment possible	
800 MV	Gas blender ²⁾	Integrated	
800 »Air« ³⁾	Flowmeter unit ¹⁾ ; operation with third gas, namely air; can be switched to mixture of O ₂ /N ₂ O or O ₂ /air	Attachment possible	For additional operating notes: see Section 11
800 V »Air« ³⁾		Integrated	

¹⁾ Setting of gas flows in L/min

²⁾ Setting of desired O₂ concentration in vol. %

³⁾ The term »Air« denotes the Romulus 800/800 V models with ancillary compressed-air unit (see Section 11 on page 19)

1 Intended Use

1.1 Applications

The Romulus 800, Romulus 800 M and Romulus 800 »Air« models are inhalation anaesthesia machines with a continuous fresh-gas flow.

The Romulus 800 V, Romulus 800 MV and Romulus 800 V »Air« models are additionally fitted with a removable automatic anaesthesia lung ventilator (Ventilog).

All models are compact, mobile machines designed for inhalation anaesthesia; they can be used in operating theatres, induction rooms and wake-up rooms.

1.2 Notes on safety

Because of the danger of explosion, never grease or oil valves on oxygen cylinders and pressure reducers for oxygen or touch them with greasy fingers.

Oxygen cylinders must not be stored together with readily flammable materials. Filled steel cylinders must not be exposed to direct sunlight or stored in the immediate vicinity of a radiator.

Valve handwheels must only be turned by hand. Never use tools!

The cylinder valves are precision parts which can easily be damaged if undue

force is applied. Valves which leak and/or do not move easily must be repaired in a workmanlike manner.

Knurled connections are only intended for manual operation.




Fig. 1 Romulus 800

39 651

Key to Figs. 1 and 2

- 1 Trolley
- 2a Cabinet 8 H
- 2b Cabinet 4 H
- 3 Hose holder
- 4 Instrument housing
- 5 Retaining arm for Pulmomat 19 or 19.K1
- 6¹⁾ Vapor 19.1
- 7a Depositing tray 1 B
- 7b Depositing tray 0.5 B
- 8 Sphygmomanometer/anaesthesia timer combination
- 9¹⁾ Oxycom® 100 D O₂ meter and monitor
- 9.1 O₂ sensor for Oxycom 100 D
- 10¹⁾ Minute Volumeter® 3000
- 11¹⁾ Precom® (airway pressure gauge with alarm)
- 12¹⁾ Circle system 7a/8 ISO
- 13¹⁾ Secretion aspirator, comprising
- 13.1 - drive (in illustration: vacuum drive)
- 13.2 - cylinder battery





Fig. 2 Romulus 800 V with integrated Ventilog

39 628

- 14¹⁾ Ventilog® anaesthesia lung ventilator
- 14.1 Pneumatic switching valve for Ventilog
- 15¹⁾ Barolog® A (airway pressure monitor)
- 15.1 Barolog measurement connection
- 16 Base plate 0.5 B

¹⁾ Device or accessory with own operating manual




39 640

Fig. 3 Romulus 800 M with blender; Ventilog additionally latched on

Key to Figs. 3 and 4
17 Blender

- 18 Writing surface
- 19 Catheter basket
- 20 Bag for manual ventilation
- 21¹⁾ Spirolog® 1 N (tidal volume monitor)
- 22¹⁾ Microbe filter (in inspiration branch)

¹⁾ Device or accessory with own operating manual



39 639

Fig. 4 Romulus 800 MV with blender and integrated Ventilog




Fig. 5 Functional schematic of Romulus 800/800 V

2 Design and Function

The numbers indicated in this section refer only to the functional schematics in Figs. 5 and 6.

All Romulus 800 models are operated with oxygen and nitrous oxide. Supply can either be effected from a central gas supply unit (referred to in the following as CS) or from gas cylinders. In the case of CS operation the shut-off valves 3 must be opened. When using gas cylinders, the cylinder pressure is indicated on the pressure gauges 1 and reduced in the pressure reducers 2 to 5 bar. The check valves 4 prevent overflow from

one gas cylinder to another or from the cylinders into the CS.

The oxygen pressure (see Fig. 5 for Romulus 800/800 V) is monitored by the O₂ deficiency signal 6 which gives an audible alarm if the O₂ pressure drops to below 2 bar. Should the pressure continue to fall, the N₂O supply is interrupted by the N₂O cutoff 7.

With the Romulus 800 M/800 MV models (Fig. 6) the N₂O cutoff is integrated into the blender 7a. The O₂ and N₂O pressure is monitored by the O₂/N₂O deficiency signal 6a; if the O₂ or N₂O pres-

sure drops below 2 bar, an audible alarm is given. The indicator 18 shows which gas has failed.

The pressure reducers 5 (Romulus 800/800 V, Fig. 5) reduce the oxygen and nitrous-oxide pressure to 1.5 bar.

The flow control valves 8 make it possible to meter the two gas flows. The flow rates can be read off in each case from two series-connected flowmeters 9 of the flowmeter unit. The two gases are routed together and passed to the Vapor selector switch 17.

Alarm and safety Devices

3.1 O₂ deficiency signal, N₂O cutoff

All Romulus 800 models are fitted with an O₂ deficiency signal and an N₂O cutoff. The 800 M/800 MV models have an additional N₂O deficiency signal.

The O₂ deficiency signal (with the 800 M/800 MV models also the N₂O deficiency signal) is designed such that an audible alarm, which cannot be reset, is given if a minimum supply pressure is dropped below.

The possible gas supply statuses are listed for all models (with the exception of Romulus 800 »Air«/800 V »Air« which are discussed in Section 11.2) in Table 1.

Explanatory notes on Table 1

A Normal operation

Oxygen and nitrous oxide are available at the prescribed pressure (see Technical Data). The O₂ deficiency signal and the nitrous-oxide cutoff are ready for operation.

Should the O₂ supply pressure drop below 2 bar, the audible O₂ deficiency alarm sounds for at least 7 seconds¹⁾. If the O₂ pressure drops below roughly 1.4 bar, the N₂O supply is reduced – until cutoff is effected – such that the pre-selected O₂ concentration is not dropped below (flow reduction). The N₂O supply is blocked completely at an O₂ pressure of ≤ 0.4 bar.

¹⁾ The anaesthetic apparatus must have been connected for at least 20s to a supply pressure ≥ 2.7 bar.

O₂ and N₂O metering can be effected again when the O₂ pressure in the system has increased to at least 2.7 bar; at this level the N₂O cutoff is also ready for further operation. The O₂ deficiency signal is ready for operation again when the O₂ pressure in the system has reached at least 2.7 bar.

O₂ failure with 800 M/800 MV models

In contrast to the 800/800 V models, the supply of N₂O is cut off completely when the alarm sounds. The green indicator for the O₂ operating pressure goes out. O₂ and N₂O metering can be effected again when O₂ supply pressure has increased to at least 2.7 bar; at this level the N₂O cutoff is also ready for further operation.

In the event of N₂O failure, O₂ can still be metered and the O₂ deficiency signal is ready for operation. No audible alarm is given. N₂O metering can be effected again when the N₂O supply pressure has reached at least 2.7 bar.

N₂O failure with 800 M/800 MV models

In contrast to the 800/800 V models, an audible alarm is given in the event of N₂O failure as well. The green indicator for the N₂O operating pressure goes out.

Moreover, should the N₂O pressure drop below 2 bar, the blender switches to 100% O₂ (slight change in flow); this switch is only made after the alarm has sounded. N₂O metering can be effected again when the N₂O supply pressure has reached at least 2.7 bar.

D O₂ and N₂O failure

Should both gases fail, the devices essentially react as described under B.

When a malfunction occurs and/or the supply pressure fluctuates outside the prescribed values, operation of the de-

Case	Model	O ₂ supply	N ₂ O supply	O ₂ deficiency signal N ₂ O deficiency signal ¹⁾ ¹⁾ only 800 M/800 MV	N ₂ O cutoff	Device operable
A	All	1	1	Ready for operation	Ready for operation	Ready for operation
B	800/800 V	0	1	Audible alarm (O ₂ ≤ 2 bar)	Reduction in N ₂ O supply (O ₂ ≤ 1.4 bar), N ₂ O cutoff (O ₂ ≤ 0.4 bar)	No, O ₂ failure
	800 M/800 MV	0	1	Audible alarm Additional visual indication (O ₂ ≤ 2 bar)	N ₂ O cutoff (O ₂ ≤ 2 bar)	
C	800/800 V	1	0	No alarm	Ready for operation No effect	No, N ₂ O failure
	800 M/800 MV	1	0	Audible alarm Additional visual indication (N ₂ O ≤ 2 bar)	N ₂ O cutoff Automatic switch to 100% O ₂	
D	800/800 V	0	0	Audible alarm	N ₂ O cutoff (O ₂ ≤ 0.4 bar)	No, O ₂ and N ₂ O failure
	800 M/800 MV	0	0	Audible alarm Additional visual indication for O ₂ and N ₂ O	N ₂ O cutoff (O ₂ ≤ 2 bar)	

Table 1 Behaviour of Romulus 800/800 V/800 M/800 MV in the event of gas failure

0 $\hat{=}$ not adequate

1 $\hat{=}$ adequate

ble and only re-continued when a fully functional compressed-gas supply has been re-established (see Section 6.5).

3.2

Other safety features

The flow control valves are provided with a guard to prevent unintentional

and shape of their rotary knobs are such that they are clearly assigned to the respective gases.

All Romulus 800 models are (optionally) provided with an O₂ flush. Turning the self-resetting lever causes an O₂ flow of roughly 55 L/min (depending on O₂ supply pressure) to be added to the fresh gas flow.

and can be used in areas where there is an explosion hazard. On electrically conductive floors the devices are protected against static charges.

They are provided with an earthing pin for potential equalization which, on non-conductive floors and in special applications, enables the apparatus to be incorporated into the room's potential equalization system (see Section 4.8).

4 Initial Preparation

4.1 Gas supply

The Romulus 800 models can be supplied with gas in the following ways:

- Connection of 3 litre cylinders (Figs. 7-10)
Attach cylinders 24 and 25 to corresponding connectors 26 and 27 and tighten connections by means of spanner 23.
- Connection of 11 litre cylinders (Figs. 8 and 9)
Insert cylinders 29 and 30 into cylinder holders 28 and clamp in position. In accordance with Figs. 8 and 9 screw high-pressure spiral tubes 31 and 32 onto corresponding cylinder valves and connectors 33 and 34 at anaesthetic apparatus and tighten with spanner (different threads for O₂ and N₂O!).




Fig. 7 Romulus 800/800 M - Back view; fitted with 2 small cylinders

- 23 Spanner
- 24 O₂ cylinder, 3 litres
- 25 N₂O cylinder, 3 litres
- 26 O₂ connector (3 litre cylinder)
- 27 N₂O connector (3 litre cylinder)




Fig. 8 Romulus 800/800 M - Back view; fitted with 2 small and 2 large cylinders

- 28 Cylinder holder for 11 litres cylinders
- 29 O₂ cylinder, 11 litres
- 30 N₂O cylinder, 11 litres
- 31 O₂ high-pressure spiral tube
- 32 N₂O high-pressure spiral tube
- 33 O₂ connector (cylinder supply, 11 litre cylinder)
- 34 N₂O connector (cylinder supply, 11 litre cylinder)




Fig. 9 Romulus 800 V/800 MV – Back view; fitted with 2 small and 2 large cylinders; Ventilog supplied with compressed air from central supply unit (CS)




Fig. 10 Romulus 800 V/800 MV – Back view; fitted with 2 small cylinders; Ventilog supplied with O₂ from anaesthetic apparatus

- 35 O₂ connecting hose (from CS)
- 36 N₂O connecting hose (from CS)
- 37 O₂ connector
- 38 N₂O connector

- 39 O₂ shut-off valve
- 40 N₂O shut-off valve
- 41 O₂/air connection at Ventilog
- 42 O₂/air connecting hose (from CS) for Ventilog

- 43 O₂ connection hose (to Ventilog)
- 44 O₂ plug-in coupling (for O₂ supply of Ventilog from anaesthetic apparatus)

– **Connection to central supply unit (CS)** (Figs. 8 and 9).

First connect connecting hoses 35 and 36 to corresponding connectors 37 and 38 (Fig. 9); then press CS plugs into appropriate wall outlet valves. Ensure that shut-off valves 39 and 40 are opened for CS operation (open by turning in anti-clockwise direction). To avoid leakage losses, the valves 39 and 40 are to be kept closed when the cylinders are open and the CS connection is not being used.

Caution:

Even when operating the apparatus from a central supply unit the cylinders should remain in position as standby supply. In this way, it is only necessary to open the cylinder valves in order to switch rapidly to cylinder supply in the event of CS failure. Check valves prevent the backflow of gas out of the cylinders into the CS.

4.2 Ventilog

4.2.1 Gas supply

If used, the Ventilog can be driven with compressed air or oxygen. The use of

oxygen is permitted if there is no compressed-air facility.

Operation with compressed air (Fig. 9)

The compressed air is to be taken from the central supply unit. The appropriate connecting hose 42 is to be screwed on to the connector 41 („Air/O₂“) on the back of the Ventilog and the plug is to be connected to the central supply unit.

In the event of compressed-air failure, no alarm is given.

Operation with O₂ (Figs. 7 and 10)

Screw connection hose 43 on to back of Ventilog („Air/O₂“ connection) and insert other end of hose (with plug connection) into coupling 44 on anaesthetic apparatus.

In the event of O₂ failure the O₂ deficiency signal in the anaesthetic apparatus sounds.

4.2.2. Pneumatic switching valve (Fig. 11)

Attach pneumatic switching valve 14.1 to circle system 7a/8 ISO (Fig. 11 or 12). Connect control line 45 to connection 46 of Ventilog and attach to pneumatic switching valve (on back). Connect cor-

rugated hoses 47 and 47a in accordance with Fig. 11. The reservoir bag 48 can also be directly attached to the pneumatic switching valve. For further information consult „Ventilog“ operating manual.

4.3 Circle system 7a/8 ISO (Figs 11 and 12)

The circle system is attached to the mounting lug of the hinged arm 49 and locked in position; its height can be adjusted. Screw on mixed-gas hose 51 at fresh-gas outlet 50 of anaesthetic apparatus and mixed-gas inlet 52 of circle system.

Pay attention to instructions given in operating manuals for circle system and ancillary equipment used, such as „Minute Volumeter 3000“, „Precom“, „Barolog A“, „Spirolog 1 N“ and „Oxycom 100 D“.

Note

The DGA¹⁾ Recommendations specify the use of an oxygen meter (e. g. Oxycom 100 D) in the inspiration section.

¹⁾ DGA = Deutsche Gesellschaft für Anästhesie und Intensivmedizin (German Association for Anaesthesiology and Intensive Care)

Secretion aspirator

(Fig. 12)

The secretion aspirator is equipped for either vacuum or ejector operation.

The following applies when employing vacuum operation:

Connect appropriate connecting hose (with connector) to aspirator and insert connector into vacuum outlet valve of central supply unit.

The following procedure is to be adopted when using an ejector-driven aspirator:

Attach bacterial filter 58 (Fig. 12) to drive-gas outlet of ejector.

Position secretion jar set with jars 60 and 61 on mount provided on base of trolley. Attach connection hose 57 to thin nozzle of jar cap 59 and other end to socket at vacuum or ejector drive; attach aspiration hose 54 to thicker nozzle of jar cap and fit secretion sight glass 53 to end of hose. The secretion sight glass can be held in position on the hinged arm 49 by means of a clamp.

For flushing the catheter, the empty jar 61 is to be filled with rinsing liquid.

For further information consult operating manual »Secretion aspirator for anaesthetic apparatus«.




Fig. 11 Hose connections between anaesthetic apparatus and Ventilog as well as between circle system and pneumatic switching valve (schematic, see also Fig. 12).

Key to Figs. 11 and 12

12	Circle system 7a/8 ISO	52	Fresh-gas inlet at circle system
14	Ventilog	53	Secretion sight glass
14.1	Pneumatic switching valve	54	Secretion aspiration hose
45	Control line	55	Shut-off valve
46	Connection for control line	56	Vent valve
47	Corrugated hose	57	Vacuum connection hose
47a	Corrugated hose	58	Bacterial filter
48	Reservoir bag	59	Jar cap with overflow safeguard and relief valve
49	Hinged arm	60	Secretion jar
50	Fresh-gas outlet at anaesthetic apparatus	61	Rinsing jar
51	Fresh-gas hose		

4.5

Microbe filter

The microbe filter 22 (Fig. 4) can be inserted between the inspiration valve and inspiration hose of the circle system using a connecting sleeve. It can also be fitted between the expiration valve and expiration hose.

Pay attention to the information given in the operating manual »Microbe filter«.




Fig. 12 Romulus 800 MV with circle system 7a/8 ISO and secretion aspirator

Vapor 19.1 (Fig. 13)

Prior to attaching the Vapor 19.1, a check must be made to ensure that the sealing rings 65 have been fitted and are in perfect condition. Mount Vapor on base 66 with plug-in adaptor 63. Locking lever 62 must face forwards. Gas leakage at the plug connection is prevented by the sealing rings 65 being pressed together by the weight of the Vapor 19.1. After mounting the Vapor 19.1, the locking lever 62 must be moved to the left until it engages, in order to ensure secure attachment to the anaesthetic apparatus.

If no Vapor is attached, the valves in the socket pins 64 form a leak-proof seal and connect the flowmeter unit and fresh-gas outlet so that O₂/N₂O gas mixtures can be metered without inhalation anaesthetic as well.

Caution

During transportation, as well as when mounting and removing the Vapor 19.1, care is to be taken to ensure that the maximum permissible angle of tilt of 45° is not exceeded.

The Vapors can be changed at any time (including during anaesthesia), since the gas flow is maintained in the metering system even when the Vapor is removed.




Fig. 13 Vapor 19.1 (plug-in system)

- | | | | |
|----|-----------------|----|--|
| 62 | Locking lever | 66 | Base (plug-in system) |
| 63 | Plug-in adaptor | 67 | Support |
| 64 | Socket pin | 68 | Whistle for O ₂ deficiency signal |
| 65 | Sealing ring | 69 | Vapor selector switch |

4.7

Anaesthetic gas removal

(Fig. 14)

The circle system and Ventilog are each equipped with an exhaust gas socket 72 or 73 to which either anaesthetic filter equipment (corrugated hose with anaesthetic filter) or an anaesthetic gas scavenging system (anaesthetic exhaust hoses 70a, b with Y-piece 71, Fig. 14) can be connected. If a Ventilog is not being used, the anaesthetic filter equipment or anaesthetic exhaust hose 70a is to be connected only to the exhaust gas socket of the circle system.

Pay attention to the instructions given in the operating manual "Anaesthetic gas scavenging system".




Fig. 14 Anaesthetic gas scavenging system on circle system and Ventilog

Key to Fig. 14

- 70a Anaesthetic gas scavenging hose
- 70b Anaesthetic gas scavenging hoses
- 71 Y-piece
- 72 Exhaust gas socket at circle system
- 73 Exhaust gas socket at Ventilog

4.8

Potential equalization

If potential equalization is required (e. g. on non-conductive floors or when monitors are latched on), it is to be provided by connecting up the cable 8301349 between the contact on the trolley of the anaesthetic apparatus and the appropriate room contact.

4.9

Bag for manual ventilation

In accordance with DGAI Recommendations¹⁾ a bag for manual ventilation 20 (Fig. 4) is to be suspended from the anaesthetic apparatus.

Recommendation:

Dräger-Bag Resutator or Dräger-Laerdal Resu Bag (suitable for sterilization in an autoclave).

¹⁾ see footnote on page 10

Testing Readiness for Operation

After cleaning, disinfection and sterilization the apparatus is always to be checked for completeness and proper operation.

A Check List for checking the apparatus prior to every start-up is given in Section 12 (information on how to use this Check List is given on page 22).

5.1 Gas supply

Check all connections for supply of O₂ and N₂O (from cylinders and from CS) on back of apparatus for tightness and freedom from leaks.

Ensure that all supply hoses are in perfect working order (visual inspection).

5.1.1

O₂ supply

(Fig. 15 for Romulus 800/800 V).

(Fig. 16 for Romulus 800 M/800 MV)

- Slowly open valve of large O₂ cylinder (11 litres). Check O₂ supply on pressure gauge 79¹⁾.
- Close cylinder valve and reduce pressure reading on gauge by opening flow control valve 77 or 86.
- Slowly open valve of small O₂ cylinder (3 litres). Check O₂ supply on pressure gauge¹⁾.

¹⁾ The O₂ cylinders are full if the O₂ pressure gauge registers 200 bar. At this pressure 11 litre cylinders and 3 litre cylinders contain 2200 or 600 litres of gas respectively.

- With the N₂O supply shut off (cylinder and/or CS) open flow control valve 77 (for O₂) or 86 and check whether flow can be set over entire measuring range of appropriate flowmeters 75 or 85; make sure that floats move freely. For this test the setting of the blender (handwheel 84, Fig. 16) on the 800 M/800 MV model has no significance (arbitrary setting).
- After reducing pressure close O₂ cylinder and flow control valve again.
- Connect apparatus to central O₂ supply (Fig. 9). Open O₂ shut-off valve 39. With 800 M/800 MV model check whether green O₂ indicator 82 lights (Fig. 16).
- Open flow control valve 77 or 86 and check whether gas flows through the appropriate flowmeters.
- Close flow control valve.
- Leave O₂ shut-off valve 39 at CS connection open.




Fig. 15 Romulus 800 - Controls




Fig. 16 Romulus 800 M (or 800 MV) - Controls

Key to Figs. 15 and 16

- 74 Flowmeter unit (Romulus 800/800 V)
- 75 O₂ flowmeters
- 76 N₂O flowmeters
- 77 O₂ flow control valve
- 78 N₂O flow control valve
- 79 Pressure gauge (cylinder pressure O₂)
- 80 Pressure gauge (cylinder pressure N₂O)

- 81 O₂ flush (Bypass)
- 17 Gas blender O₂/N₂O (Romulus 800 M/800 MV)
- 82 O₂ operating pressure indicator
- 83 N₂O operating pressure indicator
- 84 Handwheel on blender
- 85 Flowmeters: mixed gas O₂/N₂O
- 86 Flow control valve: mixed gas O₂/N₂O

5.1.2

N₂O supply

(Fig. 15 for Romulus 800/800 V,

Fig. 16 for Romulus 800 M/800 MV)

- Slowly open valve of large N₂O cylinder (11 litres). Check N₂O supply on pressure gauge 80²⁾.
- Close cylinder valve and reduce pressure reading on the gauge by opening flow control valve 78 or 86.
- Slowly open valve of small N₂O cylinder (3 litres). Check N₂O supply on pressure gauge²⁾.

²⁾ The N₂O pressure is roughly 50 bar as long as there is liquid nitrous oxide in the cylinder. When gas is extracted and the temperature drops as a result, the cylinder pressure may from time to time drop below 50 bar. The liquid content of the nitrous-oxide cylinders can only be determined by weighing them. A full 11 litre cylinder (8 kg of nitrous oxide) or 3 litre cylinder (2.25 kg of nitrous oxide) contains roughly 4000 or 1125 litres respectively.

- Open flow control valve 78 (for N₂O) or 86 and check whether flow can be set over entire measuring range of appropriate flowmeters 76 or 85; make sure that floats move freely. With 800 M/800 MV model set blender to 30 vol. % O₂.
- After reducing pressure close N₂O cylinder and flow control valve again.
- Connect apparatus to central N₂O supply (Fig. 9). Open N₂O shut-off valve 40. With Romulus 800 M/800 MV check whether green N₂O indicator 83 lights (Fig. 16).

- Open flow control valve 78 or 86 and check whether gas flows through the appropriate flowmeters.
- Close flow control valve.
- Leave N₂O shut-off valve 40 at CS connection open.

5.2

Gas deficiency signal, Blender function

5.2.1

Models 800/800 V only (Fig. 15): O₂ deficiency signal, N₂O cutoff

- Open O₂ and N₂O supply (cylinders or CS).
- Set O₂ flow and N₂O flow to 1 L/min and 2 L/min respectively.
- Shut-off O₂ supply (cylinder valve "closed" or pull CS plug).

After a brief period the O₂ deficiency alarm must sound and continue to sound for at least 7 seconds. The N₂O flow must also drop to zero (N₂O cutoff).

- Close flow control valves. Re-establish gas supply.

Gas deficiency signal, blender function

- Open O₂ and N₂O supply (cylinders or CS).
- Set 30 vol. % O₂ on blender; set mixed-gas flow to 4 L/min.
- Shut off O₂ supply (cylinder valve »closed« or pull CS plug).

After a brief period the O₂ indicator 82 must go out and the gas deficiency alarm must sound for at least 7 seconds. The flow must drop to zero (N₂O cutoff).

- Re-establish O₂ supply; shut off N₂O supply. After a brief period the N₂O indicator 83 must go out and the gas deficiency alarm must sound. The flow must however remain virtually the same (approx. 4 L/min), since the blender automatically switches to 100 vol. % O₂. The increase in the O₂ concentration is to be checked on the Oxycom 100 D oxygen meter.
- Check blender function:
Re-establish O₂ and N₂O supply. Set mixed gas flow to 4 L/min. The O₂ concentration is to be checked at the blender settings 30 and 80 vol. % (handwheel 84). For testing purposes, unscrew sensor 9.1 (Fig. 1) of Oxycom 100 D and detach mixed-gas hose from connection 52 (Fig. 11). Route gas flow directly to sensor. The O₂ indication on the Oxycom 100 D must be 30 ± 5 vol. % or 80 ± 12 vol. % O₂.
- Screw sensor and mixed-gas hose back on again.

5.3

O₂ flush (Bypass)

When the O₂ flush lever 81 (Fig. 15) is actuated, a steady flow of gas out of the

released, lever must return to initial position.

5.4

Circle system 7a/8 ISO (Fig. 11)

Check that mixed-gas hose 51 is properly connected; check that hose is in perfect condition (visual inspection). Perform functional check of circle system in accordance with appropriate operating manual.

5.5

Ventilog

- Check connection between Ventilog (connection 41, Figs. 9 and 10) and compressed-gas supply (from CS or plug-in coupling 44 of anaesthetic apparatus).
- In accordance with Fig. 11, check all connections between Ventilog, pneumatic switching valve, circle system and reservoir bag.
- Perform functional check of Ventilog as per appropriate operating manual.

5.6

Secretion aspirator

- With vacuum-driven secretion aspirator: make connection to central supply unit.
- Perform functional check of secretion aspirator in accordance with operating manual »Secretion aspirator for anaesthetic apparatus«.

5.7

Microbe filter (Fig. 4)

Check condition and installation of filter 22 in accordance with recommendation given in operating manual »Microbe filter«.

- Check level and top up if necessary.
- Perform functional check of Vapor in accordance with appropriate manual.

5.9

Anaesthetic gas removal

(Fig. 14)

Check whether anaesthetic gas scavenging unit 70a, 70b is connected to circle system (socket 72) and – if applicable – to socket 73 of Ventilog. Insert connector (at hose 70a) into coupling of central scavenging unit: this starts up the unit and the indicator at the scavenging coupling must be green.

If use is not being made of an anaesthetic gas scavenging unit, an anaesthetic filter with corrugated hose must be connected to the appropriate sockets on the circle system and – if applicable – the Ventilog. Check whether filters have been renewed. The filters must be firmly seated in the rubber collar.

5.10

Bag for manual ventilation

In accordance with DGA¹¹ Recommendations the anaesthetic apparatus must be provided with a bag for manual ventilation 20 (Fig. 4). Check function of the bag for manual ventilation by pumping manually: when the bag is squeezed, there must be an audible and perceptible stream of air out of the mask taper. When released, the bag must rapidly re-assume its original shape. If the mask taper is sealed (e. g. with a finger), it must only be possible to squeeze the bag slightly by hand.

¹¹ see footnote on page 10

Operational Use

6.1

Possible ventilation modes

Prior to connection of the circle system to the patient, the desired gas flow and O₂/N₂O mixing ratio are to be set at the flow control valves or blender. The following modes are possible: **spontaneous breathing, manual ventilation** and – if use is made of a Ventilog – **automatic ventilation**.

Spontaneous breathing

Set lever of switching valve 89 (Fig. 17) at circle system 7a/8 ISO such that it points vertically downwards (»spontaneous breathing«). The patient can thus exhale freely via check valve 90. The relief valve 88 has no function.

The selector switch 87 on the Ventilog must be in the »man/spont.« setting. Sufficient filling of the reservoir bag 48 is to be guaranteed by way of an adequate mixed-gas supply.

Manual ventilation

Set lever of switching valve 89 such that it points vertically upwards (»manual ventilation«). The air exhaled by the patient can escape via the relief valve 88 and the check valve. The airway pressure must be set at the relief valve 88. The selector switch 87 on the Ventilog must be in the »man/spont.« setting. Ventilation is effected manually by way of the breathing bag 48 with care being taken to ensure that the bag is adequately filled.

Automatic ventilation

Set lever of switching valve 89 such that it is horizontal (»automatic ventilation«). The relief valve and check valve have no function. Inspiration and expiration are effected only via the Ventilog. The selector switch 87 of the Ventilog must be in the »Autom.« setting. The Ventilog is to be set in accordance with the appropriate operating manual.

For further information consult the operating manuals »Ventilog« and »Circle System 7a/8 ISO«.

6.2

Secretion aspirator (Fig. 17)

Attach aspiration catheter to secretion sight glass 53. Open shut-off valve 55, set vacuum at vent valve 56 and extract secretion. Following aspiration, suck rinsing liquid through system. Close shut-off valve. Secretion jar must be emptied at the latest when the 600 ml mark is reached. Overflow is prevented by means of an overflow safeguard.




Fig. 17 View of circle system, Ventilog and secretion aspirator

- 87 Selector switch (Ventilog)
- 88 Relief valve
- 89 Switching valve (circle system)
- 90 Check valve

Pay attention to instructions given in operating manual »Secretion Aspirator for Anaesthetic Apparatus«.

6.3

Vapor 19.1

All Romulus 800 models can be fitted with two Vapor 19.1 vaporizers of which one can be selected in each case using the selector switch 69 (Fig. 13).

The Vapor 19.1 is switched on by pressing the locking button »zero« on the handwheel; the desired concentration is set by turning the handwheel.

To remove, turn locking lever 62 approximately 100° to the right. Lift off Vapor vertically and remove. Put replacement Vapor in position in reverse order.

For further information consult operating manual »Vapor 19.1«.

6.4

O₂ flush (Bypass)

Actuation of the O₂ flush 81 causes roughly 55 L/min of O₂ to flow into the circle system bypassing the flowmeters and Vapor.

Caution!

Observe pressure in breathing system! If the O₂ flush is actuated in an uncontrolled manner, impermissibly high airway pressures may result particularly in conjunction with automatic ventilation.

6.5

O₂ deficiency signal

When the O₂ deficiency alarm sounds, the O₂ supply is to be immediately re-established:

- With CS operation it is often sufficient to re-insert the CS plug. In the event of CS failure, the O₂ cylinder valve is to be opened and the O₂ shut-off valve 39 closed.
- In the event of O₂ deficiency during cylinder operation either the small O₂ standby cylinder is to be opened or a rapid cylinder change must be effected.

Note:

As an emergency gas supply, we recommend procuring the Dräger cylinder battery O₂/N₂O (for ordering data, see Prospectus 5303e).

N₂O deficiency

An N₂O deficiency does not lead to an audible alarm with the 800/800 V models, but it can be seen from the flowmeters. The situation is to be remedied in a manner similar to that described in Section 6.5. With the 800 M/800 MV models an alarm is given in the event of N₂O deficiency as well.

Bag for manual ventilation

The bag for manual ventilation 20 (Fig. 4) suspended from the anaesthetic apparatus in accordance with DGA¹⁾ Recommendations is intended for emergency ventilation. Information on this topic is given in the operating manuals »Bag Resutator« or »Resu Bag«.

¹⁾ See footnote on page 10.

7 Shut-Down Actions

Close O₂ and N₂O shut-off valves 39, 40 (Fig. 9). If the patient connection is detached, the following components must also be shut down:

- Switch off Vapor 19.1 (handwheel in »zero« setting);
- Close flow control valves 77 and 78 (or 86 with 800 M/800 MV models);
- Close shut-off valve 55 of secretion aspirator;
- Set switch on Ventilog to »man/spont.«.

In the event of lengthy periods of non-use (e. g. overnight), the cylinder valves must be closed and the plugs of the CS hoses detached from the wall outlet valves.

To prevent the gas deficiency alarm from sounding when shutdown is effected, the system is to be depressurized by opening the flow control valves until the flowmeter floats have dropped to their lowest position:

- With 800/800 V models, first open N₂O flow control valve, then O₂ flow control valve

- With 800 M/800 MV models first interrupt N₂O supply, open flow control valve and wait for N₂O deficiency alarm; only then is the O₂ supply to be interrupted. The O₂ pressure is also to be reduced by opening the flow control valve.
- Finally, the flow control valves are to be closed.

8 Care

8.1 Stripping down

- Detach all plug connections at CS outlet valves. The CS hoses can remain in position on the apparatus, but the CS plugs must not lie on the floor (in particular during disinfection in the Aseptor).
- Unscrew mixed-gas hose 51 from circle system.
- If anaesthetic apparatus is equipped with a Ventilog, then the following work is to be performed:
Detach all hoses between Ventilog and circle system at Ventilog (see Fig. 11). Unscrew the Ventilog switching valve from circle system.
- Remove anaesthetic gas scavenging unit or anaesthetic filter equipment from exhaust gas sockets.
- Remove circle system.
- Remove ancillary equipment which is not to be disinfected in the Dräger Aseptor (pay attention to instructions given in pertinent manuals).

- Detach vacuum connection hose of aspirator from drive and remove secretion jar set.
- Remove secretion jar (with secretion aspiration hose) and rinsing jar from holder and empty.
- Remove patient system from Ventilog (see operating manual »Ventilog«).
- The Vapors remain in position on the apparatus.

8.2 Hygiene

The upkeep of ancillary equipment must be carried out in accordance with the specifications given in the respective operating manuals. This applies for example to the circle system, Vapor 19.1, Ventilog and secretion jar set.

As regards the anaesthetic apparatus itself, the following measures apply with respect to upkeep:

8.2.1 Cleaning

Dirt on the anaesthetic apparatus is to be removed using a damp cloth soaked in standard detergent (wetting agent). Care is also to be taken to ensure that drawers, pull-out writing surfaces and the like are not forgotten. When cleaning the flowmeter unit, pay attention to the enclosed compatibility list for Plexiglas.

After cleaning, the anaesthetic apparatus and component parts are to be allowed to dry. The cylinder jackets and caps are to be removed and cleaned if there is a suspicion of corrosion and in any case at least every six months. When doing so, the condition of the gas cylinders is to be checked (visual inspection).

8.2.2 Disinfection

Disinfection in Dräger Aseptor®

Disinfection of the Romulus 800 models is to be carried out in accordance with

Aseptor«. Prior to disinfection, the Vapors must be in the zero setting (hand-wheel on »zero«). Disinfection measures for ancillary equipment, the Ventilog and the circle system are given in the respective operating manuals.

Wiping or spraying with liquid disinfectant

Wiping or spraying should only be employed if there is no possibility of using a Dräger Aseptor for disinfection purposes. Such methods only serve to reduce the number of bacteria on the surface of the device and can thus not be unreservedly recommended.

Spraying may also lead to controls sticking (switches, rotary knobs etc).

Please observe the enclosed compatibility

the flowmeter unit.

Disinfection in Dräger Purfactor®

Breathing hoses, reservoir bag, circle system, tube, mask, secretion jar, jar cap (not including ping-pong ball), aspiration hoses etc are placed on the Purfactor washing frame. All anaesthesia materials are cleaned, disinfected and dried in »Programme I«. Thermally instable materials, such as those made of PVC, are treated in »Programme II«.

8.2.3 Sterilization

Only the following items are suitable for sterilization in steam up to 120°C:

The parts of the circle system and Ventilog which carry breathing air, the secre-

tion jar) and the aspiration hoses.

Maximum temperature for secretion sight glass: 134°C.

The permissible sterilization temperature for secretion jars and rinsing jars is printed on them.

Caution!

Sterilization in steam accelerates the natural aging of rubber parts. They must thus be checked for leaks at frequent intervals.

8.3 Assembly

The disassembled parts are assembled in reserve order of stripping down (see Section 8.1). A functional check as per Section 5 »Testing Readiness for Operation« is then to be performed.

9 Servicing

To ensure that the anaesthetic apparatus is always ready for use and fully functional, we recommend concluding a servicing agreement with the Technical Customer Service of Drägerwerk AG, thus guaranteeing regular checking

and the necessary adjustments and spare part replacement. The apparatus should be serviced twice a year by the Technical Customer Service of Drägerwerk AG.

In this respect please note the informa-

tion given under »Important Notice« page 2.

The limit switch in the blenders of the Romulus 800 M/800 MV must be replaced by the Technical Customer Service of Drägerwerk AG every 4 years.

10 Technical Data

The data marked with *) refer (additionally) to the Romulus 800 »Air«/800 V »Air« special versions as described in Section 11.

Ambient temperature
for operation of Romulus: 15...35°C

Gas supply

- From central supply unit (CS) with following requirements:

Pressure O₂: 2.7 ... 5.5 bar
N₂O: 2.7 ... 5.5 bar
*) Air: 2.7 ... 5.5 bar

Flow rates

O₂ max. 20 L/min for O₂ metering
max. 24 L/min for ejector
approx. 55 L/min for O₂ flush
at 5 bar
approx. 35 L/min for O₂ flush
at 2.7 bar
approx. 30 L/min for Ventilog
(80⁺¹⁰ L/min peak flow)

N₂O: max 15 L/min for N₂O
metering

*) Air: max 15 L/min
(without Ventilog)

Screw connections on apparatus end
for O₂: M 12 x 1 mm female thread
for N₂O: M 14 x 1 mm female thread
*) for »Air«: M 20 x 1.5 mm male
thread

O₂ and N₂O connections with shut-off
valves.

- Additionally, cylinder supply option
with following possibilities:

- One 11 litre cylinder each for O₂
and N₂O
- One 3 litre cylinder each for O₂ and
N₂O
- One 11 litre cylinder each and one
3 litre cylinder each for O₂ and N₂O

High-pressure reducer integrated into
apparatus. Gas purity requirements
in accordance with European
Pharmacopoeia.

- Plug-in coupling (self-closing) for
driving Ventilog with oxygen from anaesthetic
apparatus; output pressure
same as O₂ supply pressure.

Gas metering units Models 800/800 V and

*) Models 800 »Air«/800 V »Air«

- Gas metering with flow control valves
Colour code: blue¹⁾ for O₂
grey¹⁾ for N₂O
yellow¹⁾ for »Air«

¹⁾ for Federal Republic of Germany

Different knurling of adjustment
knobs
Adjustment range approx. 5 turns

- *) Selector switch »Air/N₂O«
(manually operated)

- 2 O₂ flowmeters (connected in series)
Measuring range: 0.1 - 2 L/min
2.5 - 15 L/min

- 2 N₂O flowmeters (connected in series)
Measuring range: 0.05 - 1 L/min
1.25 - 10 L/min

- *) 1 Air flowmeter
Measuring range: 0.8 - 15 L/min

TOTAL PRESSURE

- ± 10% of set value,
- + 15/-5% with smallest scale value of O₂ flowmeter,
- 15/+5% with smallest scale value of N₂O flowmeter,
- + 15/-5% with smallest scale value of »Air« flowmeter

Gas metering units

Models 800 M/800 MV

- Gas mixing with O₂/N₂O blender
 - Adjustment range 30-100 vol. % O₂
 - Mixing accuracy ± 10% of set value
 - Inherent consumption 1.5 L/min O₂ (discharged on back)
- Gas metering with 1 flow control valve; adjustment range approx. 5 turns
- 2 Mixed-gas flowmeters (connected in series)
 - Measuring range: 0.5 - 2 L/min
 - 2.5 - 20 L/min
- Measurement accuracy (at 20°C and ± 10% of set value, + 15/-5% with smallest measuring range value and a mixture of 30% O₂ and 70% N₂O

Anaesthetic metering

2 Vapor 19.1 vaporizers for halothane or enflurane with plug-in system and selector switch. The connections of the plug-in system are automatically bypassed and sealed off with respect to the atmosphere when the Vapor is removed.

- Vapor 19.1 for halothane:
 - Adjustment range 0.2 - 4 vol. %
- Vapor 19.1 for enflurane:
 - Adjustment range 0.2 - 5 vol. %

For technical data of Vapor 19.1 see appropriate operating manual.

Gas deficiency signal/N₂O cutoff Models 800/800 V

+ Models 800 »Air«/800 V »Air«

If an O₂ supply pressure of 2.1 ± 0.1 bar is dropped below, an audible alarm, which cannot be disconnected, is given for at least 7 seconds. Should the O₂ pressure continue to drop to below roughly 1.4 bar, the N₂O supply is throttled and blocked completely at an O₂ pressure of approximately 0.4 bar.

+) In contrast to the 800/800 V models, the air supply to the »Air« flowmeter is automatically released with the 800 »Air«/800 V »Air« models should the O₂ pressure drop below 1 bar (in the system).

When the O₂ supply pressure starts to increase again, the alarm is ready for further operation as of 2.7 bar.

Models 800 M/800 MV

Should the O₂ or N₂O supply pressure drop below 2.1 ± 0.1 bar, an audible alarm is given in each case. Indicators show which gas has failed. In the event of an O₂ deficiency, the supply of N₂O is automatically cut off. Should the N₂O supply fail, the blender switches to 100 vol. % O₂ once the deficiency signal has sounded.

O₂ flush (bypass)

Approx. 55 L/min O₂ at 5 bar O₂ supply pressure

Approx. 35 L/min O₂ at 2.7 bar O₂ supply pressure

Self-resetting, no increase in pressure at Vapor.

Circle system

- For example modified circle system such as Dräger Circle System 7a or 8 ISO.
- For equipment and technical data see operating manual »Circle System 7a/8 ISO«.
- Connection for mixed-gas hose:
 - Male thread M 16 x 1.5 mm

- Drive:

By means of vacuum from central supply unit (connection via male thread ISO 228 G 1/4 A) (previous designation R 1/4")

Or by means of O₂ ejector from anaesthetic apparatus.

- Vacuum setting:

By means of vent valve, 0 to approx. -0.9 bar, vacuum indication on drive by means of pressure gauge, measuring range 0 to -1 bar.

- Effective capacity of secretion and rinsing jar: 700 mL

- Overflow safeguard, relief valve

- Rapid vent valve (only for vacuum drive)

- Bacterial filter (only for O₂ drive)

- For technical data see »Secretion aspirator« operating manual.

Ventilog anaesthesia lung ventilator

- Supply by means of O₂ plug-in coupling at anaesthetic apparatus or O₂ central supply 2-6 bar or compressed-air central supply 2-6 bar.

- Drive gas must be dry and free from oil.

- Drive gas consumption 30 L/min, peak flow 80^{+10} L/min O₂ or compressed air.

- For equipment and technical data see »Ventilog« operating manual.

Dimensions	Model 800	Model 800 V
	800 M	800 MV
	800 »Air«	800 V »Air«
- Width:	700 mm	900 mm
- Height:	1150 mm	1150 mm
- Depth:	600 mm	600 mm

Weight

Not including cylinders and circle system approx. 75 kg approx. 90 kg

Special Versions

Romulus 800 »Air«/800 V »Air«

Section 11 contains additional information necessary for operation of the Romulus 800 models »Air«. Operation of these models does however presuppose knowledge of the preceding sections of this operating manual for the standard versions.

The »Air« models differ from the standard versions in that they have an additional gas supply (compressed air) from a CS and a flowmeter unit extended to include compressed air with a selector switch for »Air« and N₂O (Figs. 18a and 18b).




Fig. 18a Romulus 800 »Air« with ancillary compressed-air unit, selector switch in setting »O₂/Air«




Fig. 18b Romulus 800 »Air« with ancillary compressed-air unit, selector switch in setting »O₂/N₂O«

- 91 Selector switch »Air«/N₂O
- 92 »Air« flowmeter
- 93 »Air« flow control valve

11.1

Design and function (Fig. 19)

The numbers in this section refer only to the functional schematic illustrated in Fig. 19.

The »Air« models are driven with oxygen, nitrous oxide and compressed air. Supply is effected from a CS or – however only for O₂ and N₂O – from gas cylinders.


For CS operation the shut-off valves 3 must be opened. When effecting supply from gas cylinders, the cylinder pressure is indicated on the pressure gauges 1 and reduced to 5 bar at the pressure reducers 2. The check valves 4 prevent overflow between the cylinders and the CS.

As regards N₂O and compressed air, it is only possible to meter one gas or the other. For this purpose the selector switch 19 is to be moved to the »N₂O« or »Air« setting.

The oxygen pressure is monitored by the O₂ deficiency signal 6 which triggers an audible alarm when the O₂ pressure drops below 2 bar. Should the O₂ pressure continue to drop, the N₂O supply is interrupted by the N₂O cutoff 7 irrespective of the selector switch setting; the apparatus switches automatically (at compressed air valve 7b) to »Air«. If the selector switch 19 is in the »Air« setting, the supply of compressed air remains connected even in the event of oxygen deficiency.

The pressure reducers 5 reduce the pressure of the oxygen and nitrous oxide to 1.5 bar.

The flow control valves 8 make it possible to meter the two gas flows (O₂ and N₂O) and the flow rate can be read off in each case from two series-connected flowmeters 9 of the flowmeter unit. As an alternative to N₂O, air can be metered via an individual flowmeter. The gases are routed together and passed to the Vapor selector switch 17.



- | | | | |
|----|----------------------------------|----|--|
| 1 | Cylinder pressure gauge | 10 | Vapor plug-in system |
| 2 | Pressure reducer 200/5 bar | 11 | O ₂ flush (Bypass) |
| 3 | Shut-off valve | 12 | Aspiration ejector (optional) |
| 4 | Check valve | 13 | Vent valve |
| 5 | Pressure reducer 5/1.5 bar | 14 | Plug-in coupling for Ventilog |
| 6 | O ₂ deficiency signal | 15 | Fresh-gas outlet |
| 7 | N ₂ O cutoff | 16 | Vapor 19.1 |
| 7b | Compressed-air valve | 17 | Vapor selector switch |
| 8 | Flow control valve | 19 | Selector switch »Air«/N ₂ O |
| 9 | Flowmeters | | |

Fig. 19 Functional schematic of Romulus 800/800 V with ancillary compressed-air unit

be switched into the fresh-gas flow for adding the desired anaesthetic. Both Vapors 16 are connected via plug-in systems 10 which permit rapid replacement. Gas can flow from the flowmeter unit to the fresh-gas outlet 15 even when no Vapor is connected.

The O₂ flush 11 (bypass) enables an O₂ flow of roughly 55 L/min (depending on CS pressure) to be added to the fresh-gas without affecting the pressure ratios (ejector system). The O₂ flush lever resets automatically.

If fitted, the aspiration ejector 12 is designed to generate a vacuum for operating the secretion aspirator. The vacuum generated (max. -0.9 bar) can be reduced via the vent valve 13.

The self-closing plug-in coupling 14 is intended for driving a Ventilog anaesthesia lung ventilator with oxygen.

associated alarm functions are tabulated below (Table 2) and then described.

Switch setting »Air« (O₂ and air)

Case 1

In the switch setting »Air«, air can be metered in combination with O₂. The gas supply to the N₂O metering branch is automatically blocked. The O₂ deficiency signal is ready for operation. In the case of central supply units with a low pressure (e. g. 3 ± 0.3 bar), it is advisable, when connecting the apparatus to the CS, to keep the O₂ metering valve at the flowmeter unit closed for a brief period to enable the pressure reservoir of the O₂ deficiency signal to fill up thus ensuring that an alarm will be given for at least 7 seconds as prescribed.

Case 2

In the event of O₂ failure during operation, air can still be metered. Should the O₂ pressure drop below 2 bar, the O₂ deficiency signal sounds for at least 7 seconds.

The O₂ deficiency signal is only ready for operation again when there is an O₂ pressure of at least 2.7 bar in the system.

Case 3

In the event of air failure during operation, O₂ can still be metered. The O₂ deficiency signal remains ready for operation. No audible alarm is given. Air can only be metered again when there is an air pressure of at least 2.7 bar in the system.

Case 4

In the switch setting »N₂O«, N₂O can be metered in combination with O₂. The gas supply to the air metering branch is automatically blocked. Mixtures of N₂O and air are not possible. The O₂ deficiency signal and N₂O cutoff are ready for operation.

Case 5

Should the O₂ pressure drop below 2 bar, the O₂ deficiency signal sounds for at least 7 seconds. Should the O₂ pressure drop below roughly 1.4 bar, the N₂O supply is reduced – until it is shut off completely – in such a manner that the pre-selected O₂ concentration is not dropped below. At an O₂ pressure of ≤ 0.4 bar, the supply of N₂O is shut off completely.

The supply of air is released at an O₂ pressure below roughly 1 bar: air can be metered via the »Air« flow control valve. Irrespective of this, the selector switch remains in the »N₂O« setting.

The O₂ deficiency signal and the N₂O cutoff are only ready for operation again when there is an O₂ pressure of at least 2.7 bar in the system. At this pressure the supply of N₂O is also released again and the air supply is blocked.

Case 6

In the event of N₂O failure, O₂ can still be metered. The O₂ deficiency signal continues to be ready for operation. No audible alarm is given, N₂O metering can only be effected again when there is an N₂O pressure of at least 2.7 bar in the system again.

11.2 Alarm and safety devices

The O₂ deficiency signal is designed as for the 800/800 V models: if an O₂ supply pressure of 2 bar is dropped below, an audible alarm is triggered which cannot be reset. Should the O₂ pressure drop below roughly 1 bar, the N₂O supply is blocked with the »Air« models and the supply of air released so that air can be metered in as emergency supply.

Case	Selector switch in setting	O ₂ supply	Air supply	N ₂ O supply	O ₂ deficiency signal	N ₂ O cutoff	Device operable
1	»Air«	1	1	Automatic blocked	Ready for operation	Ready for operation No effect	Ready for operation
2		0	1	Automatic blocked	Audible alarm (O ₂ ≤ 2 bar)	Ready for operation No effect	No, O ₂ failure
3		1	0	Automatic blocked	Ready for operation	Ready for operation No effect	No Air failed No alarm
4	»N ₂ O«	1	No effect	1	Ready for operation	Ready for operation	Ready for operation
5		0	Automatic switch to air	Audible alarm (O ₂ ≤ 2 bar)	Automatic blocking of N ₂ O (O ₂ ≤ 0.4 bar)	No, O ₂ failure	
6		1	No effect	0	Ready for operation	Ready for operation No effect	No N ₂ O failure No alarm

Table 2 Behaviour of Romulus 800/800 V models with ancillary compressed-air unit in the event of gas failure

0 \triangleq not adequate
1 \triangleq adequate

Initial preparation

Gas supply

The 800 «Air»/800 V «Air» models have an additional CS connector 94 for air (Fig. 20). The «Air» models are to be connected as described in Section 4. The «Air» connecting hose 95 is to be additionally connected to the CS.

11.4

Testing readiness for operation

Following cleaning, disinfection and sterilization the apparatus is always to be checked for completeness and proper functioning.

O₂ supply

As for Romulus 800/800 V (see Section 5.1).

N₂O supply

As for Romulus 800/800 V (see Section 5.1).

Compressed-air supply (Fig. 20)

- Check tightness of CS screw connection (at connector 94) and ensure that the CS-plug is firmly attached. Move selector switch 91 to «Air» setting. Check whether the full flow range can be set with the aid of the «Air» flow control valve 93 at the air flowmeter 92.
- Check that air flowmeter float moves freely.
- Open N₂O flow control valve: there must be no N₂O flow.
- Close N₂O and air flow control valve.

O₂ deficiency signal/N₂O cutoff/ switchover to air

- Set selector switch to «N₂O» setting. Set O₂ flow to 1 L/min and N₂O flow to 2 L/min.




Fig. 20 Additional compressed-air connection on Romulus 800/800 V
94 «Air» connector
95 «Air» connecting hose (from CS)

- Open «Air» flow control valve: There must be no flow of gas (air). Close O₂ cylinder valve or pull O₂ CS-plug. After a brief period the O₂ deficiency signal must sound and continue to sound for at least 7 seconds. The N₂O flow must also drop to zero and the air flow must increase. Re-establish O₂ gas supply. Compressed air must be automatically blocked and N₂O flow rate must again be 2 L/min.
- Set selector switch to «Air» setting. Set O₂ flow to 1 L/min and air flow to 2 L/min. Close O₂ cylinder valve or pull O₂ CS-plug. After a brief period the O₂ deficiency signal must sound and continue to sound for at least 7 seconds. The air flow must not change and the N₂O flow must remain on zero. Re-establish gas supply. Close flow control valves.

The checking of all other functions is to be performed in accordance with the description given for the Romulus 800/800 V models (see Section 5.).

11.5

Operational use

Flowmeter unit with selector switch (Figs. 18a and 18b)

The compressed air is metered at the flow control valve 93. The metered quantity is read off from the flowmeter 92. Air and O₂ can only be metered when the selector switch 91 is in the «Air» setting. In the «N₂O» setting only N₂O and O₂ can be metered.

Reference values for the O₂ concentration in O₂/air mixtures in the flow range 2–24 l/min are given in Table 3 (page 22).

	2		4		6		8		10		12		14		16		18		20		22		24			
	Luft	O ₂	Luft	O ₂	Luft	O ₂	Luft	O ₂	Luft	O ₂	Luft	O ₂	Luft	O ₂	Luft	O ₂	Luft	O ₂	Luft	O ₂	Luft	O ₂	Luft	O ₂		
O ₂ concentration in vol. %	21	2	0	4	0	6	0	8	0	10	0															
	30	1,8	0,2	3,5	0,5	5,3	0,7	7,1	0,9	8,9	1,1															
	40	1,5	0,5	3,0	1,0	4,6	1,4	6,1	1,9	7,6	2,4	9,1	2,9													
	50	1,3	0,7	2,5	1,5	3,8	2,2	5,1	2,9	6,3	3,7	7,6	4,4	8,9	5,1	10	6,0									
	60	1,0	1,0	2,0	2,0	3,0	3,0	4,0	4,0	5,0	5,0	6,1	5,9	7,1	6,9	8,1	7,9	9,1	8,9							
	70	0,8	1,2	1,5	2,5	2,3	3,7	3,0	5,0	3,8	6,2	4,6	7,4	5,3	6,7	6,1	9,9	6,8	11,2	7,6	12,4	8,3	13,7	9,1	14,9	
	80	0,5	1,5	1,0	3,0	1,5	4,5	2,0	6,0	2,5	7,5	3,0	9,0	3,5	10,5	4,0	12	4,6	13,4	5,1	14,9					
	90	0,3	1,7	0,5	3,5	0,8	5,2	1,0	7,0	1,3	8,7	1,5	10,5	1,8	11,2	2,0	14									
	100	0	2	0	4	0	6	0	8	0	10	0	12	0	14											

Table 3 Gas composition: air/O₂ (average values)

Further notes

on Operational Use, Shut-Down Actions, Care, Servicing, Technical Data, Check List, Order List and Parts List are given in the corresponding Sections 6-14.

12 Check List

In the Federal Republic of Germany, testing of the anaesthetic apparatus in accordance with the Check List on page 23 is mandatory. Please observe the recommendations or regulations in force in your country.

For users in the Federal Republic of Germany, use of this Check List is described in the following.

»The Check List on page 23 for inhalation anaesthesia apparatus must, on the basis of the corresponding operating manuals (for the anaesthetic apparatus and ancillary equipment), be brought into line with the type and configuration of the respective apparatus by means of deletions and/or additions. The Check List then contains the tests which are always to be performed on the respective model prior to start-up. The model designation and serial number of the apparatus in question are also to be entered.

The above-mentioned entries are to be transferred to the Check List (plastic) included with the anaesthetic apparatus using a waterproof felt-tip pen. The plastic Check List is to be attached to the anaesthetic apparatus by means of the bead chain.

Entries in the ACTUAL column and in the space provided for the date and signature are intended as an indication of performance of the respective tests. These entries are to be made in pencil and rubbed out again when the next set of tests is performed.

The plastic Check List must not be wiped over with cleaning agents and disinfectants, alcohol or similar solvents, since the entries made with a waterproof felt-tip pen are not resistant to such substances. Disinfection in the Aseptor is however permitted«.

DrägerModel _____
Serial No. _____

- Knowledge of valid operating manuals is an absolute prerequisite
- Delete where not applicable; make additions where necessary

Signature _____

Checking prior to start-up

What?	How?	Desired	Actual
Anaesthetic gas Cylinder supply	Open valves	Pressure O ₂ > 50 bar N ₂ O > 30 bar	tick off if okay
Central supply	Insert plug-in coupling Open O ₂ /N ₂ O metering valves	Indicator green Flow present
Anaesthetic gas scavenging system	Insert plug-in coupling	Indicator green
Anaesthetic filter	Condition of filter	Filter replaced
O ₂ -flush (bypass)	Actuate switch	Flow present
Vapor	Zero setting Level Selector switch	Locked Adequate Switch setting correct
Plug-in system	Connection	Plug-in system locked
Circle system	Hoses Reservoir bag Absorber Volumeter Volumeter heating Airway pressure gauge Measurement connections Valve discs (insp. and exp.) Mixed-gas hose	Completeness and tight fit
Soda lime	Condition of lime		
O ₂ meter	Functional check, calibration	Lime renewed, no colour change Functional
Monitors	Functional check, calibration	Functional
Freedom from leaks for non-rebreathing and modified circle system	Seal relief valve and Y-piece, flow 0.4 L/min	Pressure ≥ 20 mbar for 10 seconds
Relief valve	Relief valve 20 mbar, seal Y-piece, flow 10 L/min	Constant pressure 20 ± 5 mbar
System Non-rebreathing/modified circle system	Selector switch	Switch setting correct
Ventilator	Connections to circle system Switch on, check settings, seal Y-piece during inspiration	Tight Airway pressure present
Secretion aspirator	Switch on, seal aspiration hose	Vacuum present
Bag for manual ventilation, for emergency ventilation	Check completeness Check bag	Complete Functioning properly
Additions			

The articles such as masks, Y-pieces, corrugated hoses and connectors preceded by a (P) symbol correspond in terms of connection dimensions to the ISO Draft Standard ISO/DIS 5356. In addition there are a number of items which can be used both with ISO and with Dräger standard connecting elements.

Name and Description	Order Code
Basic versions	
A) Basic versions with flowmeter unit	
Romulus 800 basic version 1	M 25804 ¹⁾
For operation from two 11 litre gas cylinders and two 3 litre gas cylinders as standby or for operation from a central supply unit. With trolley, connections block with flowmeter unit, cylinder connections, O ₂ and N ₂ O high-pressure spiral tubes, hinged arm, mixed gas hose, O ₂ deficiency signal with N ₂ O cut-off	
Romulus 800 basic version 5	M 25808 ¹⁾
For operation from a central supply unit and two 3 litre gas cylinders as standby. With trolley, connections block with flowmeter unit, cylinder connections, hinged arm, mixed gas hose, O ₂ deficiency signal with N ₂ O cut-off	
B) Basic versions with gas mixer	
Romulus 800 M basic version 1	M 25731 ¹⁾
For operation from two 11 litre gas cylinders and two 3 litre gas cylinders as standby or for operation from a central supply unit. With trolley, connections block with gas blender, cylinder connections, O ₂ and N ₂ O high-pressure spiral tubes, hinged arm, mixed gas hose, O ₂ deficiency signal with N ₂ O cut-off	
Romulus 800 M basic version 5	M 25735 ¹⁾
For operation from a central supply unit and two 3 litre gas cylinders as standby. With trolley, connections block with gas blender, cylinder connections, hinged arm, mixed gas hose, O ₂ deficiency signal with N ₂ O cut-off	
C) Basic versions with flowmeter unit, envisaged for installation of a Ventillog	
Romulus 800 V basic version 3	M 25806 ¹⁾
For operation from two 11 litre gas cylinders and two 3 litre gas cylinders as standby or for operation from a central supply unit. With trolley, connections block with flowmeter unit, cylinder connections, O ₂ and N ₂ O high-pressure spiral tubes, hinged arm, mixed gas hose, O ₂ deficiency signal with N ₂ O cut-off	
Romulus 800 V basic version 7	M 25810 ¹⁾
For operation from a central supply unit and two 3 litre gas cylinders as standby. With trolley, connections block with flowmeter unit, cylinder connections, hinged arm, mixed gas hose, O ₂ deficiency signal with N ₂ O cut-off	
D) Basic versions with gas blender, envisaged for installation of a Ventillog	
Romulus 800 MV basic version 3	M 25733 ¹⁾
For operation from two 11 litre gas cylinders and two 3 litre gas cylinders as standby or for operation from a central supply unit. With trolley, connections block with gas blender, cylinder connections, O ₂ and N ₂ O high-pressure spiral tubes, hinged arm, mixed gas hose, O ₂ deficiency signal with N ₂ O cut-off	
Romulus 800 MV basic version 7	M 25737 ¹⁾
For operation from a central supply unit and two 3 litre gas cylinders as standby. With trolley, connections block with gas blender, cylinder connections, hinged arm, mixed gas hose, O ₂ deficiency signal with N ₂ O cut-off	
Anaesthesia cabinets required for all basic versions, option of:	
Cabinet 4 H	2 M 18093
Anaesthesia cabinet with one drawer and depositing tray	

¹⁾ Order codes change when the cylinder connection differs from the standard German connection. Will be inserted by Drägerwerk AG or its agency.

Cabinet 6 H	2 M 18095
Anaesthesia cabinet with three drawers and depositing tray	
Cover plates required for all basic versions, option of:	
Cover plate, with rail	M 25450
Simple plate with depositing tray	
Base plate 1 B, with rail	M 25839
(Required for latching on Dräger Monitors).	
Accessories required for operation	
Circle system 7a	M 23074
With 2 carbon dioxide absorbers, inhalation and exhalation valve, relief valve, breathing bag, 3 corrugated hoses, soda lime filler funnel, 1 set of valve discs and sealing rings	M 25690
a) For operation from cylinders:	
O ₂ cylinder, 11 litres	B 2710 ¹⁾
200 bar, contents 2200 litres of oxygen	
N ₂ O cylinder, 11 litres	B 2660 ¹⁾
Filled with 8 kg of nitrous oxide	
O ₂ cylinder, 3 litres	B 2533 ¹⁾
200 bar, contents 600 litres of oxygen	
N ₂ O cylinder, 3 litres	B 2540 ¹⁾
Filled with 2.25 kg of nitrous oxide	
Cylinder jacket, 11 litres	M 3691
With valve cap and base ring, for 11 litre cylinder	
Cylinder jacket, 3 litres	M 8035
For 3 litre cylinder	
b) For operation from a central supply unit, option of:	
O ₂ connecting hose, 3 m	M 22344
O ₂ connecting hose, 5 m	M 22345
Connection - device end:	
Clamping screw M 12 × 1	
Connection - wall end:	
Angled plug-in socket for O ₂	
N ₂ O connecting hose, 3 m	M 22350
N ₂ O connecting hose, 5 m	M 22351
Connection - device end:	
Clamping screw M 14 × 1	
Connection - wall end:	
Angled plug-in socket for N ₂ O	
Required for basic versions 3 and 7:	
Ventillog	84 04 500
For controlled ventilation during anaesthesia	
For operation from Romulus 800 M or 800 MV:	
Connecting hose, 0.6 m	M 25050
Alternatively (for operation from a central supply unit), option of:	
O ₂ -compressed-air connecting hose, 3 m	M 22494
O ₂ -compressed-air connecting hose, 5 m	M 22495
Connection accessories required:	
Connecting hose 2/1 m	84 04 758
Pneumatic switching valve	84 04 950
	84 05 276
Alternatively:	
Connecting hose 2/1 m	84 04 758
Manual switching valve	84 05 305
	84 05 295
Accessories recommended for monitoring	
1. For continuous measurement and monitoring of O ₂ in gas inhaled:	
Oxycom 100 D oxygen meter and monitor	68 03 255
Measuring range 0-100% O ₂ at normal pressure. With sensor housing, cable and sensor capsule. Connection elements required for oxygen meter:	
a) If anaesthesia timer, sphygmomanometer or combination is being used:	
O ₂ meter holder 11	M 21483
b) For screwing directly onto threaded connection of hinged arm:	
O ₂ meter holder 10	M 21478

Name and Description	Order Code	Name and Description	Order Code
2. For continuous measurement and monitoring of airway pressure in circle system:		2. Via anaesthetic filter:	
a) Barolog A	83 02 930	Anaesthetic filter equipment 2	M 21262
With visual and audible alarm in the event of disconnection or obstruction Connection accessories required: Measurement connection	M 25638	For filtering out harmful anaesthetic vapours. With 5 anaesthetic filters. 1 set of anaesthetic filter equipment in each case is required for the circle system and Ventilog	
Alternatively:		Secretion aspirator, ejector	M 26136
b) Respiratory pressure gauge with alarm "Precom"	E 9711	Secretion aspirator, vacuum	M 26137
Gives an audible alarm if a set pressure value is not attained within 15 seconds Alternatively:		Required for operation of vacuum aspirator M 26137, option of:	
c) Respiratory pressure gauge without alarm	72 64 325	Vacuum connecting hose, 3 m	M 22353
3. For continuous measurement and monitoring of minute volume, tidal volume and frequency:		Vacuum connecting hose, 5 m	M 22354
a) Spirolog 1 N	83 02 760	Connection - device end: Cap nut R 3/4"	
Connection accessories required: Spirolog sensor housing	M 26844	Connection - wall end: Angled plug-in socket for vacuum	
S-set: sensor, set of 5	84 03 735		
Alternatively:			
b) For measurement of minute volume and tidal volume: Minute Volumeter 3000	2 M 18250		
Special accessories		Ventilog, fitted with latching elements, for latching onto base plate 1 B with basic versions 1 and 5	84 05 200
Anaesthesia timer/sphygmomanometer Combined	M 14626	For operation from Romulus 800 or 800 M:	
Anaesthesia timer	M 14692	Connecting hose 1.2 m	M 25518
Sphygmomanometer	M 14691	Alternatively (for operation from central supply unit), option of:	
Blood pressure cuff, size 3 for adults	M 13790	O₂-compressed-air connecting hose, 3 m	M 22494
Blood pressure cuff, size 2 for children	M 20139	O₂-compressed-air connecting hose, 5 m	M 22495
Blood pressure cuff, size 1 for infants	M 20140	Connecting accessories required: Connection hose 2/1.5 m	84 04 732
		Pneumatic switching valve	84 04 950
		Alternatively:	84 05 276
		Connection hose 2/1.5 m	84 04 732
		Manual switching valve	84 05 305
			84 05 295
Gas-analysis measurement connections For continuous CO ₂ and O ₂ measurement during anaesthesia	M 18074	Depositing tray 0.5 B for latching on to Ventilog or a Dräger Monitor (e. g. Barolog A or Spirolog 1 N)	2M 17680
Set of microbe filters 644 St For insertion between inhalation valve and inhalation hose in circle system. 5 filters per pack. Can be sterilized 20 times	67 27 260	Depositing tray 1 B for latching on to 2 adjacent Dräger Monitors	84 07 025
► Set of microbe filters 644 St Suitable for ISO circle system, comprising:		Dummy casing 2 H 0.5 B For height compensation, for latching on next to Dräger monitor	M 25625
Set of microbe filters 644 St	67 27 260	Holder with rail For attachment to basic versions 3 and 7 on lefthand side next to anaesthesia cabinet. Designed to accommodate holder with refuse bag or the like	M 26390
ISO-Set for microbe filters	M 26930	Base plate 0.5 B for height compensation. To latch on to a built-in Ventilog in Romulus basic versions 3 and 7	M 26965
E-Vapor 19.1/5% pin safety With safety filling system	DB 01023	Writing surface	2 M 18260
H-Vapor 19.1/4% pin safety With safety filling system	DB 01022	Holder For accommodating a Pulmomat 19	M 25205
Special accessories for safety filling system		Bypass (accessory) for basic versions 1 and 5	M 25740
Filling hose, Enflurane	M 26299	Bypass (accessory) for basic versions 3 and 7	M 25741
Filling hose, Halothane	M 26297	For expanding basic units with flowmeter unit Romulus 800/800 V to form basic units Romulus 800 »Air«/800 V »Air« (air as third gas in addition to O ₂ and N ₂ O):	
E-Vapor 19.1/5%	DB 01041	Ancillary compressed-air unit	M 26597
H-Vapor 19.1/4%	DB 01040	Connecting hoses required, option of:	
		Compressed-air connecting hose, 3 m	M 23193
		Compressed-air connecting hose, 5 m	M 23235
Possibilities for simultaneous removal of anaesthetic vapours from circle system and Ventilog		Drawer inserts for anaesthesia cabinets:	
1. Via an ejector system:		Insert 2 (2 compartments)	G 12101
a) Accessory set - waste anaesthetic gas scavenging system 15	M 26094	Insert 4 (4 compartments)	G 12102
For basic versions 3 and 7 with built-in Ventilog		Insert 6 (6 compartments)	G 12103
b) Accessory set - waste anaesthetic gas scavenging system 16	M 26095	See prospectus 5325e for Babylog N anaesthesia lung ventilator for infants	
For basic versions 1 and 5 with latched-on Ventilog		Semi-open anaesthesia system 2.1 for spirometry, connection facilities for respiratory pressure gauge, Volumeter and waste anaesthetic gas scavenging system	M 23210
Further accessories:		Infants' anaesthesia kit according to Kuhn	M 14832
Plug	G 60225		M 25634
Connecting hoses required, option of:			
Anaesthetic exhaust hose, 3 m	G 60305		
Anaesthetic exhaust hose, 5 m	G 60306		

Name and Description	Order Code
Waste anaesthetic gas exhaust connection, for Kuhn kit	M 23100
Set of Rendell-Baker masks, for Kuhn kit	►M 25838
(qty. 4) 1 mask each, sizes 0-3	M 24526
Accessory set for children's circle system	M 26702
Bag Reserator	M 11900
Resu-bag	21 09 832
Hook	M 26349
Dust cover for basic versions 1 and 5	2 M 6832
Dust cover for basic versions 3 and 7	2 M 6835
Earth cable 3.2 m	83 01 349
For sealing of cylinder connections for 3 litre standby cylinders:	
O₂ screw plug	M 6620
N₂O lock nut	M 6621
Anaesthesia accessories (special prospectus 5301.0e)	
Replacement parts for sterilization	
Jar set, for secretion aspirator	M 26355
Patient set, for Ventilog including anaesthetic exhaust socket	84 05 040
Circle system 7a	M 23074
a) Pneumatic switching valve, for Ventilog	►M 25690
Connection hose 2/1 m	84 04 950
Alternatively:	►84 05 276
Connection hose 1/1.5 m	84 04 758
b) Manual switching valve, for Ventilog	84 04 732
Connection hose 2/1 m	84 05 305
Alternatively:	►84 05 295
Connection hose 2/1.5 m	84 04 758
Respiratory pressure gauge	E 9726
Further connecting accessories, eg. Y-pieces, masks, etc.: see prospectus 5301.0e «Anaesthetic Accessories».	

No in Fig. 21	Name	Item No.
2	Extension arm	M 25 410
3	Screw	M 14 075
4	Set of T-screws (Qty. 2)	M 22 191
5	Washer	M 25 419
6	T-screw	M 19 816
7	Mixed-gas hose comprising Sealing ring (Qty. 2)	M 17 734
8	Hose	
7	Set of sealing rings (Qty. 10)	M 22 189
9	T-screw	M 24 953
10	O-ring set (Qty. 10)	U 15 314
11	Corner piece	M 23 396
12	Spanner 32/22 mm	M 12 401
13	Connecting hose comprising Sealing ring	M 25 050
14	Hose	
13	Set of sealing rings (Qty. 10)	M 23 454
15	High-pressure spiral tube O ₂ comprising High-pressure spiral tube O ₂	M 7 571
16	Sealing ring (Qty. 2)	
16	Set of sealing rings (Qty. 10)	D 20 065
16a	Set of sealing rings (Qty. 10) for pin index cylinder connection (not illustrated)	M 23 498
-	Set of T-screws (Qty. 2) for pin index cylinder connection (not illustrated)	M 23 615
18	High-pressure spiral tube N ₂ O comprising High-pressure spiral tube N ₂ O	M 7 572
19	Sealing ring (Qty. 2)	
19	Set of sealing rings (Qty. 10)	M 23 439
20	Set of screws (Qty. 10)	2M 17568
21	Hose 2/1 m	84 04 758




Fig. 21 Component parts for all Romulus models (see Parts List on page 26).
 The item numbers are not identical with the item numbers in the other Figures



DRÄGERWERK AG LÜBECK
FEDERAL REPUBLIC OF GERMANY

✉ Postfach 1339, 🏠 Moislinger Allee 53/55
D-2400 Lübeck 1
☎ (451) 882-0, 📠 26807-0, FAX (451) 882-2080